



Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that the our product

e-motion M15

including the Lithium ion batteries fulfil the requirements of

IATA chapter 2.3.2.4

and therefore may be carried as checked baggage without restriction.

Condition precedent to meet the above listed requirement is the transportation in a way that an unintentionally switching on and damage by other luggage is avoided.

The batteries must be removed for transportation and carried (protected from damage and from short circuit) in the **passenger cabin**.

The pilot in command has to be informed about the location of the batteries.

The remote control ECS of the e-motion could be transported also under IATA chapter 2.3.5.9 which allows the transportation via cabin baggage.

We strongly recommend obtaining the approval by corresponding airline company before commencement date respectively departure date.

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2018**.

Albstadt, 2017-12-20

Alber GmbH

A handwritten signature in black ink, appearing to read "H. Reiner", written over a light blue horizontal line.

Hansjörg Reiner
Managing Director